

ANTISERA - SALMONELLA

SSI Diagnostica *Salmonella* antisera can identify O and H antigens and the capsular antigen Vi. The *Salmonella* antisera are for slide agglutinations or phase inversions, respectively.

SALMONELLA SEROTYPING KITS (CE)

ART. NO.	PRODUCT	NO. OF TESTS	PACKING	AVAILABILITY
18349	<i>Salmonella</i> Sero-Quick Group kit	150	1 box	In stock
18350	<i>Salmonella</i> Sero-Quick ID kit	150	1 box	In stock

SALMONELLA ANTISERA AND MONOCLONAL ANTIBODY FOR SEROLOGICAL CONFIRMATION BY SLIDE AGGLUTINATION (CE)

ART. NO.	PRODUCT	REACTS WITH FACTOR	NO. OF TESTS	PACKING	AVAILABILITY
40211	Poly A-E + Vi	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 15, 19, 20, 27, 34, 46 + Vi	150	3 mL vial	In stock
73517	Poly A-E + Vi	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 15, 19, 20, 27, 34, 46 + Vi	50	1 mL vial	In stock
44807	Poly A-I + Vi	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 19, 20, 22, 23, 24, 25, 27, 34, 46 + Vi	150	3 mL vial	In stock
73518	Poly A-I + Vi	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 19, 20, 22, 23, 24, 25, 27, 34, 46 + Vi	50	1 mL vial	In stock
48954	Poly A-S + Vi	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 27, 28, 30, 34, 35, 38, 39, 40, 41, 46 + Vi	150	3 mL vial	In stock
73519	Poly A-S + Vi	1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 27, 28, 30, 34, 35, 38, 39, 40, 41, 46 + Vi	50	1 mL vial	In stock
81061	Poly 42-67	42, 43, 44, 45, 47, 48, 50, 51, 52, 53, 54, 55, 56, 57, 58, 59, 60, 61, 62, 63, 65, 66, 67	50	1 mL vial	In stock
40290	Poly H	a b c d i z ₁₀ E G k L r z y Z ₄ 1 z ₆ z ₂₉	150	3 mL vial	In stock
73520	Poly H	a b c d i z ₁₀ E G k L r z y Z ₄ 1 z ₆ z ₂₉	50	1 mL vial	In stock
15514	Vi	Vi (monoclonal antibody)	150	3 mL vial	In stock
73521	Vi	Vi (monoclonal antibody)	50	1 mL vial	In stock

SALMONELLA O GROUP POOL ANTISERA FOR SLIDE AGGLUTINATION (CE)

ART. NO.	PRODUCT	REACTS WITH GROUP/FACTOR	NO. OF TESTS	PACKING	AVAILABILITY
40212	OMA	A, B, D, E, L 1, 2, 3, 4, 5, 9, 10, 12, 15, 19, 21, 46	150	3 mL vial	In stock
73522	OMA	A, B, D, E, L 1, 2, 3, 4, 5, 9, 10, 12, 15, 19, 21, 46	50	1 mL vial	In stock
40213	OMB	C, F, G, H 6, 7, 8, 11, 13, 14, 20, 22, 23, 24, 25	150	3 mL vial	In stock
73523	OMB	C, F, G, H 6, 7, 8, 11, 13, 14, 20, 22, 23, 24, 25	50	1 mL vial	In stock
40214	OMC	I, J, K, M, N 16, 17, 18, 28, 30	150	3 mL vial	In stock
73524	OMC	I, J, K, M, N 16, 17, 18, 28, 30	50	1 mL vial	In stock
40215	OMD	O, P, Q, R, S, T, U, V, W 35, 38, 39, 40, 41, 42, 43, 44, 45	150	3 mL vial	In stock
73525	OMD	O, P, Q, R, S, T, U, V, W 35, 38, 39, 40, 41, 42, 43, 44, 45	50	1 mL vial	In stock
40216	OME	X, Y, Z 51, 52, 47, 48, 50, 51, 52	150	3 mL vial	In stock

ANTISERA - SALMONELLA

SALMONELLA O GROUP POOL ANTISERA FOR SLIDE AGGLUTINATION (CONT.)

ART. NO.	PRODUCT	REACTS WITH GROUP/FACTOR	NO. OF TESTS	PACKING	AVAILABILITY
73526	OME	X, Y, Z 51, 52, 47, 48, 50, 51, 52	50	1 mL vial	In stock
40217	OMF	53, 54, 55, 56, 57, 58, 59	150	3 mL vial	In stock
73527	OMF	53, 54, 55, 56, 57, 58, 59	50	1 mL vial	In stock
40218	OMG	60, 61, 62, 63, 65, 66, 67	150	3 mL vial	In stock
73528	OMG	60, 61, 62, 63, 65, 66, 67	50	1 mL vial	In stock
40219	O3 Complex	3, 10, 15, 19, 34	150	3 mL vial	In stock
73529	O3 Complex	3, 10, 15, 19, 34	50	1 mL vial	In stock

SALMONELLA O GROUP ANTISERA FOR SLIDE AGGLUTINATION

ART. NO.	PRODUCT	INCL. IN GROUP	NO. OF TESTS	PACKING	AVAILABILITY
40220	O:2	A	150	3 mL vial	In stock
23839	O:4	B	150	3 mL vial	In stock
40221	O:6, 7, 8	C ₁ , C ₂ , C ₃ , H	150	3 mL vial	In stock
65335	O:7, 8	C ₁ , C ₂ , C ₃	150	3 mL vial	In stock
40224	O:8	C ₂ , C ₃	150	3 mL vial	In stock
40276	O:9	D ₁ , D ₂ , D ₃	150	3 mL vial	In stock
40226	O:3, 10, 15	E ₁ , E ₄	150	3 mL vial	In stock
23843	O:1, 3, 19	E ₁ , E ₄	150	3 mL vial	In stock
40227	O:11	F	150	3 mL vial	In stock
40228	O:13, 22, 23	G ₁ , G ₂	150	3 mL vial	In stock
40229	O:6, 14, 24, 25	H	150	3 mL vial	In stock
40230	O:16	I	150	3 mL vial	In stock
40231	O:17	J	150	3 mL vial	In stock
40232	O:18	K	150	3 mL vial	In stock
40233	O:21	L	150	3 mL vial	In stock
40234	O:28	M	150	3 mL vial	In stock
40235	O:30	N	150	3 mL vial	In stock
40236	O:35	O	150	3 mL vial	In stock
40237	O:38	P	150	3 mL vial	In stock
40238	O:39	Q	150	3 mL vial	In stock
40239	O:40	R	150	3 mL vial	In stock
40240	O:41	S	150	3 mL vial	In stock
40241	O:42	T	150	3 mL vial	In stock
40242	O:43	U	150	3 mL vial	In stock
40243	O:44	V	150	3 mL vial	In stock
40244	O:45	W	150	3 mL vial	In stock
40245	O:47	X	150	3 mL vial	In stock
40247	O:48	Y	150	3 mL vial	In stock
40248	O:50	Z	150	3 mL vial	In stock
40249	O:51	51	150	3 mL vial	In stock
40250	O:52	52	150	3 mL vial	In stock
40251	O:53	53	150	3 mL vial	In stock
40252	O:54	54	150	3 mL vial	In stock
40253	O:55	55	150	3 mL vial	In stock
40254	O:56	56	150	3 mL vial	In stock
40255	O:57	57	150	3 mL vial	In stock
40257	O:58	58	150	3 mL vial	In stock

ANTISERA - SALMONELLA

SALMONELLA O GROUP ANTISERA FOR SLIDE AGGLUTINATION (CONT.) (CE)

ART. NO.	PRODUCT	INCL. IN GROUP	NO. OF TESTS	PACKING	AVAILABILITY
40258	O:59	59	150	3 mL vial	In stock
40259	O:60	60	150	3 mL vial	In stock
40260	O:61	61	150	3 mL vial	In stock
40261	O:62	62	150	3 mL vial	In stock
40263	O:63	63	150	3 mL vial	In stock
40264	O:65	65	150	3 mL vial	In stock
40265	O:66	66	150	3 mL vial	In stock
40266	O:67	67	150	3 mL vial	In stock
40267	T1		150	3 mL vial	In stock
40268	T2		150	3 mL vial	In stock

SALMONELLA O FACTOR ANTISERA FOR SLIDE AGGLUTINATION (CE)

ART. NO.	PRODUCT	INCL. IN GROUP	NO. OF TESTS	PACKING	AVAILABILITY
40269	O:1	A, B, D ₁ , D ₂ , D ₃ , E ₄ , G, H, R, T, V, X, 51, 53, 59	150	3 mL vial	In stock
40272	O:5	B	150	3 mL vial	In stock
40273	O:6 ₁	C ₁ (6 ₁ +), C ₂	150	3 mL vial	In stock
40274	O:7	C ₁	150	3 mL vial	In stock
40277	O:10	E ₁	150	3 mL vial	In stock
57543	O:12	A, B, D ₁ , D ₃	150	3 mL vial	In stock
40278	O:14	C ₁ , H, K	150	3 mL vial	In stock
40279	O:15	E ₁	150	3 mL vial	In stock
40280	O:19	E ₄	150	3 mL vial	In stock
40281	O:20	C ₂ , C ₃	150	3 mL vial	In stock
40282	O:22	G	150	3 mL vial	In stock
40283	O:23	G	150	3 mL vial	In stock
65334	O:24, 25	H	150	3 mL vial	In stock
40285	O:25	H	150	3 mL vial	In stock
40286	O:27	B, D ₃	150	3 mL vial	In stock
40287	O:34	E ₁	150	3 mL vial	In stock
40288	O:46	D ₂ , D ₃	150	3 mL vial	In stock

SALMONELLA H PHASE POOL ANTISERA FOR SLIDE AGGLUTINATION (CE)

ART. NO.	PRODUCT	REACTS WITH FACTOR	NO. OF TESTS	PACKING	AVAILABILITY
40291	HMA	a + b + c + d + i + z ₁₀	150	3 mL vial	In stock
73530	HMA	a + b + c + d + i + z ₁₀	50	1 mL vial	In stock
40292	HMB	E + G	150	3 mL vial	In stock
73531	HMB	E + G	50	1 mL vial	In stock
40293	HMC	k + L + r + z + y + Z ₄	150	3 mL vial	In stock
73532	HMC	k + L + r + z + y + Z ₄	50	1 mL vial	In stock
40294	HMD	1,2 + 1,5 + 1,6 + 1,7 + z ₆	150	3 mL vial	In stock
73533	HMD	1,2 + 1,5 + 1,6 + 1,7 + z ₆	50	1 mL vial	In stock
40295	HME	Z ₂₉ , Z ₃₅ , Z ₃₈ , Z ₃₉ , Z ₄₁ , Z ₄₂	150	3 mL vial	In stock
73534	HME	Z ₂₉ , Z ₃₅ , Z ₃₈ , Z ₃₉ , Z ₄₁ , Z ₄₂	50	1 mL vial	In stock
40296	HMF	Z ₂₇ , Z ₃₆ , Z ₅₂ , Z ₅₃ , Z ₅₄ , Z ₅₅ , Z ₅₆ , Z ₅₇	150	3 mL vial	In stock
73535	HMF	Z ₂₇ , Z ₃₆ , Z ₅₂ , Z ₅₃ , Z ₅₄ , Z ₅₅ , Z ₅₆ , Z ₅₇	50	1 mL vial	In stock
48953	HMG	Z ₄₄ , Z ₆₀ , Z ₆₁ , Z ₆₄ , Z ₆₅ , Z ₆₇ , Z ₆₈	150	3 mL vial	In stock
73536	HMG	Z ₄₄ , Z ₆₀ , Z ₆₁ , Z ₆₄ , Z ₆₅ , Z ₆₇ , Z ₆₈	50	1 mL vial	In stock
40297	H:L	All phases having l	150	3 mL vial	In stock
73537	H:L	All phases having l	50	1 mL vial	In stock

ANTISERA - SALMONELLA

SALMONELLA H PHASE POOL ANTISERA FOR SLIDE AGGLUTINATION (CONT.)

ART. NO.	PRODUCT	REACTS WITH FACTOR	NO. OF TESTS	PACKING	AVAILABILITY
40298	H:E	All phases having e	150	3 mL vial	In stock
73538	H:E	All phases having e	50	1 mL vial	In stock
40299	H:G	All phases having g and m,t	150	3 mL vial	In stock
73539	H:G	All phases having g and m,t	50	1 mL vial	In stock
40366	H:Z ₄	All phases having z ₄	150	3 mL vial	In stock
73540	H:Z ₄	All phases having z ₄	50	1 mL vial	In stock
22667	H:q,s,t,p,u	q, s, t, p, u	150	3 mL vial	In stock
73542	H:q,s,t,p,u	q, s, t, p, u	50	1 mL vial	In stock
91024	H:G ÷ H:m,t	All phases having g (not m,t) Used to distinguish between H:g,m,t and H:m,t	50	1 mL vial	In stock

SALMONELLA H PHASE ANTISERA FOR SLIDE AGGLUTINATION

ART. NO.	PRODUCT	NO. OF TESTS	PACKING	AVAILABILITY
40300	H:a	150	3 mL vial	In stock
23847	H:b	150	3 mL vial	In stock
40301	H:c	150	3 mL vial	In stock
23849	H:d	150	3 mL vial	In stock
53844	H:e,h	150	3 mL vial	In stock
53824	H:e,n,x	150	3 mL vial	In stock
53825	H:e,n,z ₁₅	150	3 mL vial	In stock
23850	H:i	150	3 mL vial	In stock
40303	H:k	150	3 mL vial	In stock
40304	H:r	150	3 mL vial	In stock
40305	H:y	150	3 mL vial	In stock
40306	H:z	150	3 mL vial	In stock
40307	H:z ₆	150	3 mL vial	In stock
40308	H:z ₁₀	150	3 mL vial	In stock
40309	H:z ₂₉	150	3 mL vial	In stock
40310	H:z ₃₅	150	3 mL vial	In stock
40311	H:z ₃₆	150	3 mL vial	In stock
40312	H:z ₃₈	150	3 mL vial	In stock
40313	H:z ₃₉	150	3 mL vial	In stock
40314	H:z ₄₁	150	3 mL vial	In stock
40315	H:z ₄₂	150	3 mL vial	In stock
40316	H:z ₄₄	150	3 mL vial	In stock
40317	H:z ₅₂	150	3 mL vial	In stock
40318	H:z ₅₃	150	3 mL vial	In stock
40319	H:z ₅₄	150	3 mL vial	In stock
40320	H:z ₅₅	150	3 mL vial	In stock
40321	H:z ₅₇	150	3 mL vial	In stock
40609	H:z ₆₀	150	3 mL vial	In stock
40322	H:z ₆₁	150	3 mL vial	In stock
40323	H:z ₆₄	150	3 mL vial	In stock
40324	H:z ₆₅	150	3 mL vial	In stock
40325	H:z ₆₇	150	3 mL vial	In stock
52233	H:z ₆₈	150	3 mL vial	In stock
40326	H:z ₇₁	150	3 mL vial	In stock
40327	H:z ₈₁	150	3 mL vial	In stock
40328	H:z ₈₃	150	3 mL vial	In stock
57544	H:z ₉₁	150	3 mL vial	In stock

ANTISERA - SALMONELLA

SALMONELLA H FACTOR ANTISERA FOR SLIDE AGGLUTINATION

ART. NO.	PRODUCT	NO. OF TESTS	PACKING	AVAILABILITY
40331	H:f	150	3 mL vial	In stock
40329	H:h	150	3 mL vial	In stock
40332	H:m	150	3 mL vial	In stock
40333	H:p	150	3 mL vial	In stock
40334	H:q	150	3 mL vial	In stock
40335	H:s	150	3 mL vial	In stock
40336	H:t	150	3 mL vial	In stock
40337	H:u	150	3 mL vial	In stock
40338	H:v	150	3 mL vial	In stock
40339	H:w	150	3 mL vial	In stock
40330	H:x	150	3 mL vial	In stock
40340	H:2	150	3 mL vial	In stock
40341	H:5	150	3 mL vial	In stock
40342	H:6	150	3 mL vial	In stock
40343	H:7	150	3 mL vial	In stock
40344	H:z ₁₃	150	3 mL vial	In stock
40345	H:z ₁₅	150	3 mL vial	In stock
40346	H:z ₁₆	150	3 mL vial	In stock
40347	H:z ₂₃	150	3 mL vial	In stock
40348	H:z ₂₄	150	3 mL vial	In stock
40349	H:z ₂₈	150	3 mL vial	In stock
40350	H:z ₃₂	150	3 mL vial	In stock
40351	H:z ₅₁	150	3 mL vial	In stock

SALMONELLA H:R-PHASE ANTISERA FOR SLIDE AGGLUTINATION

ART. NO.	PRODUCT	NO. OF TESTS	PACKING	AVAILABILITY
40352	H:Rj	150	3 mL vial	In stock
40353	H:Rz ₂₇	150	3 mL vial	In stock
40354	H:Rz ₃₃	150	3 mL vial	In stock
52234	H:Rz ₃₄	150	3 mL vial	In stock
40355	H:Rz ₃₇	150	3 mL vial	In stock
40356	H:Rz ₄₀	150	3 mL vial	In stock
40357	H:Rz ₄₃	150	3 mL vial	In stock
40358	H:Rz ₄₅	150	3 mL vial	In stock
40359	H:Rz ₅₆	150	3 mL vial	In stock
48961	H:Rz ₅₉	150	3 mL vial	In stock
48959	H:Rz ₆₆	150	3 mL vial	In stock
40360	H:Rz ₇₂	150	3 mL vial	In stock
40610	H:Rz ₇₄	150	3 mL vial	In stock
40361	H:Rz ₇₆	150	3 mL vial	In stock
40362	H:Rz ₇₇	150	3 mL vial	In stock
52235	H:Rz ₇₉	150	3 mL vial	In stock

ANTISERA - SALMONELLA

SALMONELLA SG POOL ANTISERA FOR PHASE INVERSION (S. GARD METHOD¹) ©

ART. NO.	PRODUCT	REACTS WITH FACTOR	NO. OF TESTS	PACKING	AVAILABILITY
40363	SG 1	a + b + c + z ₁₀	30	3 mL vial	In stock
40364	SG 2	d + i + e,h	30	3 mL vial	In stock
40365	SG 3	k + y + L*	30	3 mL vial	In stock
40367	SG 5	e,n,x + e,n,z ₁₅	30	3 mL vial	In stock
40368	SG 6	1* + z ₆	30	3 mL vial	In stock
40369	SGG	All phases having g or m,t	30	3 mL vial	In stock
40370	SGZ4	z ₄₁ z ₂₃ + z ₄₁ z ₂₄ + z ₄₁ z ₃₂	30	3 mL vial	In stock

¹Gard,S. (1938) Z. Hyg.: 120, 615 L*: All phases having l

1*: 1,2,1,5,1,6,1,7

SALMONELLA H PHASE ANTISERA FOR PHASE INVERSION (S. GARD METHOD¹) ©

ART. NO.	PRODUCT	NO. OF TESTS	PACKING	AVAILABILITY
86865	H:i	30	3 mL vial	In stock
40381	H:r	30	3 mL vial	In stock
40382	H:z	30	3 mL vial	In stock
40383	H:z ₂₉	30	3 mL vial	In stock
53950	H:Rz ₃₄	30	3 mL vial	In stock
40384	H:z ₃₅	30	3 mL vial	In stock
40385	H:z ₃₆	30	3 mL vial	In stock
40386	H:z ₃₈	30	3 mL vial	In stock
40387	H:z ₃₉	30	3 mL vial	In stock
40388	H:z ₄₁	30	3 mL vial	In stock
40389	H:z ₄₂	30	3 mL vial	In stock
40390	H:z ₄₄	30	3 mL vial	In stock
40391	H:z ₅₂	30	3 mL vial	In stock
40392	H:z ₅₃	30	3 mL vial	In stock
40393	H:z ₅₄	30	3 mL vial	In stock
40394	H:z ₅₅	30	3 mL vial	In stock
40395	H:z ₅₇	30	3 mL vial	In stock
48957	H:Rz ₅₉	30	3 mL vial	In stock
40396	H:z ₆₀	30	3 mL vial	In stock
40397	H:z ₆₁	30	3 mL vial	In stock
40398	H:z ₆₄	30	3 mL vial	In stock
40399	H:z ₆₅	30	3 mL vial	In stock
48955	H:Rz ₆₆	30	3 mL vial	In stock
48958	H:z ₆₇	30	3 mL vial	In stock
53951	H:z ₆₈	30	3 mL vial	In stock
40401	H:z ₈₁	30	3 mL vial	In stock
40402	H:z ₈₃	30	3 mL vial	In stock

¹Gard,S. (1938) Z. Hyg.: 120, 615

